



NAME 200901-20200831-HNMR  
EXPNO 1  
PROCNO 1  
Date\_ 20200831  
Time 14.32  
INSTRUM spect  
PROBHD 5 mm PABBO BB/  
PULPROG zg30  
TD 32768  
SOLVENT CDCl3  
NS 8  
DS 0  
SWH 8012.820 Hz  
FIDRES 0.244532 Hz  
AQ 2.0447731 sec  
RG 203  
DW 62.400 usec  
DE 6.50 usec  
TE 298.0 K  
D1 1.00000000 sec  
TD0 1

===== CHANNEL f1 =====  
SFO1 400.1424710 MHz  
NUC1 1H  
P1 13.35 usec  
SI 65536  
SF 400.1400014 MHz  
WDW EM  
SSB 0  
LB 0.30 Hz  
GB 0  
PC 1.00